

Musterstudie 2009

Subject Information

Subject Name: Mike MM Mustermann	Subject ID: 54321
Gender: Male	Age: 36
Height: 179 cm	Regression Set: ERS93/Polgar
Weight: 76 kg	BMI: 23,7

Test Information

Test Date: 02.06.2009 17:39	Values at BTPS
Device ID: 0000111106	Device Software Ref: 40913 V1.01
Number of Blows: 3	Number of Good Blows: 3
Best FEV1 within: 0,51 L	Best FEV6 within: 0,54 L

Results

Parameter	Predicted	Test 1	Test 2	Test 3	Best	% Pred
FEV1 (L)	4,16	4,25	3,74	3,45	4,25	102
FEV6 (L)	5,03	4,64	4,10	3,79	4,64	92
FEV1/FEV6 (ratio)	0,81	0,92	0,91	0,91	0,92	114

Obstructive Index



COPD (GOLD) Classification

FEV1 Meas./FEV1 Pred. \geq 80%

FEV1 Meas./FEV1 Pred. $<$ 80%

FEV1 Meas./FEV1 Pred. $<$ 50%

FEV1 Meas./FEV1 Pred. $<$ 30%



FEV1/FEV6 \geq 0.70

FEV1/FEV6 $<$ 0.70 & FEV1 \geq 80% Pred.

FEV1/FEV6 $<$ 0.70 & FEV1 $<$ 80% Pred.

FEV1/FEV6 $<$ 0.70 & FEV1 $<$ 50% Pred.

FEV1/FEV6 $<$ 0.70 & FEV1 $<$ 30% Pred.


Interpretation



Lung Age: 30

Obstructive Index: 102%

Interpretation: Normal, not COPD



Comments

Signature: _____

Date: _____