



the difference is clear

MDI Training Spacer

Ensure optimal and effective MDI delivery

- Enables patients to practice technique
- Allows for assessment of MDI & Spacer use
- One-way valve
- Low-cost disposable
- Fits most placebo devices

The most common delivery device for treatment of asthma is a metered dose inhaler (MDI). It is estimated that over 70% of patients do not use the MDI properly. A spacer is commonly prescribed along with the inhaler to help eliminate improper MDI use and enhance delivery of the medication to the lungs.

The low-cost, disposable, Vitalograph MDI training spacer allows the patient to practice technique and for their health care professional to assess the patients MDI and spacer use before they obtain the real thing. The MDI training spacer fits most placebo devices and has a one-way valve preventing the patient from breathing into the MDI placebo.

Vitalograph is the world's leading provider of outstanding respiratory diagnostic, screening and monitoring devices. With a heritage of excellence spanning half a century, Vitalograph continue to make pioneering contributions to effective respiratory care and enhanced quality of life.





the difference is clear



Ordering Info:

Catalogue No.

29772 MDI Training Spacer (box of 150)

Vitalograph is the world's leading provider of outstanding respiratory diagnostic, screening and monitoring devices. With a heritage of excellence spanning half a century, Vitalograph offers a complete range of pulmonary function testing products used to diagnose and manage COPD, asthma and other lung diseases.

For more information on the complete line of Vitalograph products, please visit our website at www.vitalograph.com or call our customer care team at 1-800-255-6626

Vitalograph Inc.
13310 West 99th Street
Lenexa
Kansas 66215 USA
Tel.: (913) 888-4221
Fax: (913) 888-4259
E-mail: vitcs@vitalograph.com

Vitalograph GmbH.
Rellinger Straße 64a
20257 Hamburg
Germany
Tel.: +49 40 547391-0
Fax: +49 40 547391-40
E-mail: info@vitalograph.de

Vitalograph Ltd.
Maids Moreton
Buckingham
MK18 1SW England
Tel.: +44 (0) 1280 827110
Fax: +44 (0) 1280 823302
E-mail: sales@vitalograph.co.uk

Vitalograph (Ireland) Ltd.
Gort Road Business Park
Ennis
Co. Clare Ireland
Tel.: (065) 6864100
Fax: (065) 6829289
E-mail: sales@vitalograph.ie

